

# Densah® Bur Usage Tracking Log



	VPLTT	VPLTT-S	VT1525	VT1525-S	VT1828	VT1828-S	VS2228	VT2535	VT2535-S	VT2838	VT2838-S	VS3238	VT3545	VT3545-S	VT3848	VT3848-S	VS4248	VT4555	VT4555-S	VT4858	VT4858-S	VS5258	VT5565	VT5868
1																								
2																								
3																								
4																								
5																								
6																								
7																								
8																								
9																								
10																								
11																								
12																								
13																								
14																								
15																								
16																								
17																								
18																								
19																								
20																								



## **Densah® Burs Maintenance & Replacement**

**Surgical drills and burs should be replaced after  
12-20 osteotomies\* or sooner when they are dulled or worn.**

**Cutting and Densifying must be done under constant water irrigation. A pumping motion is required to prevent overheating. It is recommended to keep a spare set of Densah® Burs on hand in the event replacement is needed during a surgery.**

\*Chacon GE, Bower DL, Larson PR, et al. Heat production by three implant drill systems after repeated drilling and sterilization. J Oral Maxillofac Surg. 2006;64(2):265-269.